

# SONICLIFE



Soniclife is the definitive keyboard against all pain from arthrosis and inflammation, even chronic in hands, wrists, elbows, neck, cervical and fibromyalgia where a recent study published in the journal on novel physiotherapies, confirms a reduction of over 75% in ultralaser is applied to the hands.

Thanks to its technologies of ultrasound, magnetotherapy, electrostimulation and ultralaser, usable with eight comfortable automatic programs, Soniclife will relieve your pain in a simple way and directly at your home.

## DUPUYTREN CONTRACTURES

De Dupuytren's disease is a disease of the hand that causes contracture of the fingers. Dupuytren's disease consists of a localized thickening, tending to extend, of the palm aponeurosis, a fascial structure placed under the skin of the palm of the hand.

## DE QUERVAIN SYNDROME

De Quervain's tendonitis is an inflammatory disease that affects the Abductor long and Extensor short tendons of the thumb as they pass through the 1st extensor canal.

## TRIGGER FINGER

The trigger finger is a very common pathology that presents with pain and evident jerk of one or more fingers of the hand in flexion and extension.

## RHIZOARTROSIS

The trapeziometacarpal joint, located at the base of the thumb, is very important because thanks to its great freedom of movement, it allows the thumb to have a very specialized action compared to the fingers, allowing opposition and consequently grip.

## CARPAL TUNNEL

Carpal tunnel syndrome is a very common disease, due to compression of the nerve inside the carpal canal.

## RADIOCARPIC ARTHROSIS

The term arthrosis defines a chronic developmental disease that affects the joints, the anatomical basis of which is represented by a degenerative process of the cartilages that cover the joint bone heads and which undergo wear phenomena.

# SONICLIFE PROGRAMS/PATHOLOGY

01 TRIGGER FINGER - 02 OSTEOARTHRISIS - 03 DUPUYTREN CONTRACTURES 04 DE QUERVAIN SYNDROME - 05 CARPAL TUNNEL SYNDROME - 06 RHIZOARTROSIS - 07 ALGODYSTROPHY - 08 FIBROMYALGIA - 9A ULTRASOUNDS - 9B ULTRASOUNDS WITH OPTIONALS - 10 MAGNETO-LIGHT



## SONICLIFE USE OF THE APPLIANCE

Remove the device and place it on a stable horizontal rigid surface, and near a power outlet. Check that the voltage of the power supply corresponds to the voltage and frequency on the appliance plate. Before inserting the power cable into the socket, make sure that the power switch is set to 0 (zero). When turned on, bringing the switch to position 1 turns on the DISPLAY on which the software presentation message appears.

At the end of the previous display, SELECT A PROGRAM appears (if you do not choose a program, the last run is proposed again)

At this point, by using the keys ▲ ▼ one of the available programs can be selected. With the same keys, the TENS is adjusted when the program has started. After choosing the program, press the TIME / START key to start the treatment, to pause it, press the same button, while keeping the TIME / START button pressed, the treatment is definitively terminated. The BEEP emits 3 notes close together, the STOP LED lights up and the DISPLAY shows the message PROGRAM ENDED. You can choose another program below.

### TECHNICAL PARAMETERS.

Power supply: 230VAC 50Hz - Power: 50W - Working time: min 10 max 45 min On the back and on the ultrasound head are the identification plates with the technical data. CE marking.



# HSD

STRADA ACQUASALATA N. 11-INT. 6 - 47899 - SERRAVALLE (RSM)  
[www.hsd.sm](http://www.hsd.sm) [info@hsd.sm](mailto:info@hsd.sm)

Tel. 0549.941496

 **NUMERO VERDE**  
**8000 200 55**